

FINDING THE GENERATOR

- Press the Therapy Off button (O) and move the sleep remote near your generator. When the sleep remote produces two beeps with the second beep at a higher tone than the first, you have found the generator location.
- If the sleep remote produces three rapid beeps, you have not found the generator location. Try again.

AUDIO TONES

TWO BEEPS WITH THE SECOND BEEP AT A HIGHER TONE THAN THE FIRST



- Communication with the generator was successful.

THREE RAPID BEEPS



- Communication with the generator was not successful, try again.
- Stimulation is at the upper or lower limit.

STATUS INDICATORS



Replace sleep remote batteries



If the light is on and not flashing, your generator battery may be low, call your doctor



If this light is on and flashing communication with your generator has been interrupted.



Stimulation strength is at the upper limit



Stimulation strength gauge

Inspire Sleep Remote™ Quick Guide



Inspire Sleep Remote Model 2500



For more information, visit www.InspireSleep.com/Patients



STATUS RING

- White**.....Therapy is off
- Green—Solid**.....Therapy is on and stimulation is on
- Green—Dimming & Brightening**.....Therapy is on, stimulation is paused

TURNING THERAPY ON

- Press the Therapy On (▶||) button and hold the sleep remote over your generator for 10 seconds or until the sleep remote produces an audio tone.
- When the therapy is turned on, the generator will produce a brief stimulation pulse.

TURNING THERAPY OFF

- Press the Therapy Off (○) button and hold the sleep remote over your generator for 10 seconds or until the sleep remote produces an audio tone.

DECREASING STIMULATION STRENGTH

- Press the Decrease button (−) and hold the sleep remote over your generator for 10 seconds or until the sleep remote produces an audio tone.
- If the sleep remote produces three beeps immediately when you press the Decrease button (−), stimulation strength is at the lower limit.

INCREASING STIMULATION STRENGTH

- Press the Increase button (+) and hold the sleep remote over your generator for 10 seconds or until the sleep remote produces an audio tone.
- If the sleep remote produces three beeps immediately when you press the Increase button (+), stimulation strength is at the upper limit.

PAUSING THERAPY

- While the therapy is on, press the Therapy On button (▶||) and hold the sleep remote over your generator for 10 seconds or until the sleep remote produces an audio tone.