#### **Finding the Generator**

- Press the Therapy Off button (O) and move the sleep remote near your generator. When the sleep remote produces two beeps, with the second beep at a higher tone than the first, you have found the generator location.
- If the sleep remote produces three rapid beeps, you have not found the generator location. Try again.

#### **Audio Tones**

11	Two beeps—the second beep at a higher tone than the first	Communication with the generator was successful
	Three rapid beeps	<ul> <li>Communication with the generator was not successful, try again</li> <li>Stimulation is at the upper or lower limit</li> </ul>
33	Two repeating beeps of the same tone	<ul> <li>Sleep remote is updating the generator</li> <li>Wait until complete &amp; the success tone is played</li> <li>Update may take several seconds</li> </ul>

#### **Status Indicators**



Manufactured For: Inspire Medical Systems, Inc. 5500 Wayzata Blvd., Suite 1600 Golden Valley, MN 55416 USA 1-844-672-4357 www.inspiresleep.com

© 2023 Inspire Medical Systems, Inc. All Rights Reserved.

200-422-102 REV B 2024-01

# Welcome to a better night's sleep.

# **Inspire Sleep Remote Quick Guide**





Inspire Sleep Remote 2580 Rx Only



Status Ring		
White — Solid	Therapy is off	
Green — Solid	Therapy is on and stimulation is active	
Green — Slowly dimming & brightening	Therapy is on and stimulation is paused or therapy is in start delay	

# **Turning Therapy ON**

- Press the Therapy On (►II) button and hold the sleep remote over your generator for 10 seconds or until you hear an audio tone.
- When therapy is turned on, the generator will produce a brief stimulation pulse.

# **Turning Therapy Off**

 $\circ$  Press the Therapy Off (  ${\rm O}$  ) button and hold the sleep remote over your generator for 10 seconds or until you hear an audio tone.



# **Decreasing Stimulation Strength**

- Press the decrease button ( ) and hold the sleep remote over your generator for 10 seconds or until you hear an audio tone.
- If you hear three beeps immediately when you press the decrease button (-), stimulation strength is at the lower limit.

#### **Increasing Stimulation Strength**

- Press the increase button (+) and hold the sleep remote over your generator for 10 seconds or until you hear an audio tone.
- If you hear three beeps immediately when you press the increase button (+), stimulation strength is at the upper limit.

#### **Pausing Therapy**

• While the therapy is on, press the Therapy On button (▶II) and hold the sleep remote over your generator for 10 seconds or until the sleep remote produces an audio tone.

# **Pairing with Inspire App**

 Remove battery cover and hold down pairing button for 3 seconds until Bluetooth light starts blinking. Follow instructions on Inspire App to complete pairing.